

Call: +48 89 527 11 30

# **DOD2 – DOG OVULATION DETECTOR**

effective mating has never been so easy





## **DOD2 – DOG OVULATION DETECTOR**



# Operating mode

#### **USB** interface

• remote software update and data transfer to a computer

#### Standard / advanced

• select operation mode

#### **Easy to disinfect**

· washed under running water, which makes it easier to keep clean

#### THE KIT INCLUDES

- · a carrying case with a foam insert,
- computer,
- a lead.



- · manual,
- 4 AA-type batteries 1.5 V (LR6),
- USB cable for communication with a





#### Technical data

Approximate device weight	380 g
Dimensions	33.0 x 8.0 x 7.5 cm
Probe length	12.5 cm
Power supply	4 AA-type batteries 1.5 V (LR6)
Battery charge indicator	graphical
Battery low indicator	automatic
Current input	11 mA to 54 mA (depending on the backlight intensity)
Measurement control	single-chip microcomputer
Approximate continuous operation time on a set of alkaline batteries	150 hours
Display	2.4" LCD with LED backlight
Membrane	keyboard
Data transmission	via USB
Update	via USB
Data record	internal memory
Memory capacity	100 animals / 200,000 measurements with date and time
Measuring range	0 to 2000 units
<b>Display resolution</b>	10 units
Recommended operating temperature	0° to 50°C
Recommended storage temperature	5° to 45°C

## MAIN CHARACTERISTICS

- · reliable results,
- · large display with energy-saving LED backlighting,
- safe for animals and humans,
- does not require periodic calibration,
- memory for 100 dogs and 200,000 measurements with date and time,
- 2-year warranty,
- data transmission to a computer for easy interpretation of results (oestrus graphs, printouts, reports, archiving).

### ADDITIONAL FEATURES

- real-time clock,
- · pop-up menu,
- · setting the display operating parameters,
- · language selection,
- · automatic shutdown time change,
- · measurement settings.

The manufacturer's website provides downloadable **free software** for firmware updates, data transmission and cataloguing on a computer.

